


InxPROFIdrive V2.1

Release Notes	
Status	RELEASED OBSOLETE
Date	2012-12-04
Summary	Several optimizations and bugfixes on profile implementation

Changes in this version
Changes in v2.1 since v2.0 ----- <ul style="list-style-type: none">- Optimization: Replace PROFIdrive Dynamic memory allocation by Static Allocation- Optimization: serialize handling of Parameter requests in case of multi-Parameter request. Just one parameter per application cycle will be processed in order to minimize the load within application cycle.- fixed Type of Parameter High-/Low-limits variables. The type of limits variables must depend on the type of Parameter value.- fixed #9: SendMailbox does not get deleted when shutting down profile/application- fixed #11: wrong order in parameter 974 initial value. Handling for parameter 974 reworked.- fixed #10: Access with subindex to non-indexed parameter- Refuse change request access to Parameter description/text- Added some functions to support PNU971, PNU977 (load and store parameter Set)

Restrictions
TBD

Known Issues					
Type	Key	Summary	FixVersion	Status	Resolution
<div style="border: 1px solid orange; padding: 10px;"> JIRA project doesn't exist or you don't have permission to view it. View these issues in JIRA</div>					
Reihenfolge: type, key, summary, fixversion, status, resolution					

Resources
■ InxPROFIdrive_v2.1.zip

[Return to overview page...](#)